

General Purpose Cannula  
(MC-14200 and 2,783 others)

## Product composition and precautions

Writing	Review	Approval

- **Components**

Needle cap, Needle, Hub

- **Precautions**

- A. This product is disposable and must be discarded after use. Re-sterilization and reuse are prohibited.
- B. Sterilization cannot be guaranteed if the individual packaging is opened or damaged, so do not use.
- C. Being pricked by a contaminated product can lead to various infectious diseases.
- D. Do not use needles that are damaged or bent prior to use.
- E. Do not use if foreign substances are found in the product prior to use.
- F. Do not use for purposes other than general procedures.