

Nanum Syringe  
(NS-50 and 50 others)

# Product composition and precautions

| Writing | Review | Approval |
|---------|--------|----------|
|         |        |          |

- **Components**

Cylinder, Plunger, Gasket, Clip (50cc only), Front cap, Discharge unit, Plunger cap (Applies only to detachable types)

- **Precautions**

1) Warning

: This product is for single use only. Dispose of after use. Do not re-sterilize or reuse.

2) General Precautions

: This product is a medical device. Be sure to familiarize yourself with the intended use, directions for use, and precautions before use.

3) Precautions for Use

- ① This product is a single-use medical device. Do not reuse.
- ② Do not use after the expiration date.
- ③ Do not use if the packaging is damaged or foreign matter is found before use.
- ④ Open immediately before use and use immediately after opening.